

## 11p Subtelomere FISH Probe (HRAS)

### Introduction

The Chromosome 11p Subtelomere FISH Probe is designed to detect presence, integrity and copy number of 11p subtelomeric chromosomal regions.

### Intended Use


To measure chromosome 11 copy number and presence of the intact human 11p subtelomeric chromosomal regions.

### Cont.

### Color

11p Subtelomere (HRAS) FISH Probe      CytoGreen

### Probe Design



STP 11p FISH Probe covers a subtelomeric chromosomal region within locus 11p15.5.

Not to Scale

### Cat. No.

### Volume

CT-STP11P-10-G

10 Tests (100 µL)

### Signal Pattern Interpretation

#### Normal Patterns

2G

#### Abnormal Patterns

Other Patterns

- 1) O'Connor C. *Nature Education*. 1(1):171 (2008).
- 2) Tsuchiya KD. *Clin Lab Med*. 31(4):525-42 (2011).
- 3) Ried T, et al. *Hum Mol Genet*. 7(10):1619-26 (1998).
- 4) Knight SJ & Flint J. *J Med Genet*. 37(6):401-9 (2000).
- 5) Jalal SM, et al. *Genet Med*. 5(1):28-34 (2003).

\* CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.

DCN032

© CytoTest Inc.

[www.cytotest.com](http://www.cytotest.com)

V2016.06.01  
T-07-10-STP11P-G-EN



CytoTest Inc.  
1395 Piccard Drive, Suite 308  
Rockville, MD 20850, USA