

17p Subtelomere FISH Probe (VPS53)

Introduction

The Chromosome 17p Subtelomere FISH Probe is designed to detect presence, integrity and copy number of 17p subtelomeric chromosomal regions.

Intended Use


To measure chromosome 17 copy number and presence of the intact human 17p subtelomeric chromosomal regions.

Cont.

Color

17p Subtelomere (VPS53) FISH Probe CytoGreen

Probe Design



STP 17p FISH Probe covers a subtelomeric chromosomal region within locus 17p13.3.

Not to Scale

Cat. No.

Volume

CT-STP17P-10-G

10 Tests (100 µL)

Signal Pattern Interpretation

Normal Patterns

2G

Abnormal Patterns

Other Patterns

- 1) O'Connor C. *Nature Education*. 1(1):171 (2008).
- 2) Tsuchiya KD. *Clin Lab Med*. 31(4):525-42 (2011).
- 3) Ried T, et al. *Hum Mol Genet*. 7(10):1619-26 (1998).
- 4) Knight SJ & Flint J. *J Med Genet*. 37(6):401-9 (2000).
- 5) Jalal SM, et al. *Genet Med*. 5(1):28-34 (2003).

* CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.

DCN032

© CytoTest Inc.



CytoTest Inc.
1395 Piccard Drive, Suite 308
Rockville, MD 20850, USA